

9300 series

Commode chairs

Our 9300 series commode chairs are truly **multi-functional**. Big wheels in front for sharper turning or big rear wheels to shift its gravity point or 4 small wheels to reduce the overall length. A reclinable backrest, adjustable legrests and removable and collapsible armrests... All the comfort you need at your fingertips.



model 9301



model 9302



model 9303



model 9300

Images for illustration purposes only. Certain features may be optional.



Opening up new horizons in mobility

Technical data 9300

Total width	620 / 660 / 690 / 730 / 760 mm
Total length	1100 mm
Total height	1070 mm
Seat width	390 / 430 / 460 / 500 / 530 mm
Height armrests	260 mm
Seat height	500 mm

Height backrest	580 mm
Seat depth	450 mm
Backrest angle	92° - 149°
Weight	25,56 kg
Users weight	100 kg
Additional info	CE

Technical data 9301

Total width	620 / 660 / 690 / 730 / 760 mm
Total length	890 mm
Total height	1070 mm
Seat width	390 / 430 / 460 / 500 / 530 mm
Height armrests	260 mm
Seat height	500 mm

Height backrest	580 mm
Seat depth	450 mm
Backrest angle	92° - 149°
Weight	26,5 kg
Users weight	100 kg
Additional info	CE

Technical data 9302

Total width	580 / 620 / 650 / 690 / 720 mm
Total length	840 mm
Total height	1070 mm
Seat width	390 / 430 / 460 / 500 / 530 mm
Height armrests	260 mm
Seat height	500 mm

Height backrest	580 mm
Seat depth	450 mm
Backrest angle	92° - 149°
Weight	24,5 kg
Users weight	100 kg
Additional info	CE

Technical data 9303

Total width	520 / 560 / 590 / 630 / 660 mm
Total length	850 mm
Total height	1070 mm
Seat width	390 / 430 / 460 / 500 / 530 mm
Height armrests	260 mm
Seat height	470 mm

Height backrest	580 mm
Seat depth	450 mm
Backrest angle	92° - 149°
Weight	23,5 kg
Users weight	100 kg
Additional info	CE

Measurements in mm, kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible

Options



L06 - Headrest



L55 - Headrest



B22 - Spreader

Your guarantee of quality

Regarding quality, the **9300 series** comply with the highest norms, laid down by certified test institutes in most European countries.

VERMEIREN India Rehab Private Ltd.

400, Red Sanders, Boulevard, Sri City (SEZ),
Satyavedu (M), Chittoor District,
Andhra Pradesh - 517646, India

Tel : +91 8576 666600

www.vermeiren.in
info@vermeiren.in

For more information, contact your local dealer



VERMEIREN

Opening up new horizons in mobility



Subject to technical alterations