


Eagle - 620 / 625

Patientlift up to 150 / 175 kg

The patient lifts Eagle - model 620 and model 625 combine functionality, safety and reliability.

The high-quality smooth-running castors in conjunction with the narrow overall width allow easy transfer even in home environments. For storage, the Eagle - model 625 can easily be folded without tools and stored upright.



Images for illustration purposes only. Certain features may be optional.



Opening up new horizons in mobility

Technical data 620 - 625

model 620

| | |
|-----------------------------|---------------|
| Total length | 1115 mm |
| Total width closed / opened | 550 / 1115 mm |
| Highest point | 1740 mm |
| Weight | 37 kg |
| Users weight | 150 kg |
| Additional information | CE |



Optional slings

VERMEIREN offers different slings for our stable and easy-to-use patient lifts. Whether as a standard version, hygiene or comfort version with upper body support; our slings offer the desired support for each persons needs.

Your guarantee of quality

Regarding quality, the patientlift **Eagle - model 620 and 625** complies with the highest norms, laid down by certified test institutes in most European countries.



model 625


| | |
|-----------------------------|---------------|
| Total length | 1290 mm |
| Total width closed / opened | 680 / 1225 mm |
| Highest point | 1840 mm |
| Weight | 42 kg |
| Users weight | 175 kg |
| Additional information | CE |

Measurements in mm, kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible

HYGIENIC SLING

available in
M - L - XL


| | A | B | C | D |
|---|-----|-----|-----|-----|
| M | 740 | 220 | 580 | 410 |
| L | 780 | 260 | 660 | 470 |
| X | 860 | 310 | 710 | 500 |



COMFORT SLING

available in
M - L - XL


| | A | B | C | D |
|---|-----|-----|-----|-----|
| M | 750 | 350 | 470 | 270 |
| L | 820 | 380 | 490 | 400 |
| X | 900 | 430 | 590 | 400 |



STANDARD SLING

available in
M - L - XL

| | A | B | C | D |
|---|-----|-----|-----|-----|
| M | 740 | 370 | 430 | 400 |
| L | 790 | 420 | 480 | 430 |
| X | 850 | 480 | 520 | 450 |



Measurements in mm, kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible

VERMEIREN India Rehab Private Ltd.

400, Red Sanders, Boulevard, Sri City (SEZ),
Satyavedu (M), Chittoor District,
Andhra Pradesh - 517646, India

Tel : +91 8576 666600

www.vermeiren.in
info@vermeiren.in

For more information, contact your local dealer



VERMEIREN

Opening up new horizons in mobility

Subject to technical alterations