

Sedna

Always mobile



The Sedna scooter is the perfect transport scooter. Lightweight, compact in size and foldable under a minute super easy to carry-on like a travelling suitcase on a boat, train or airplane (including its battery!). Ideal for short distances and handy to negotiate sharp turns in small rooms, restaurants or shops.



Images for illustration purposes only. Certain features may be optional.

Opening up new horizons in mobility

Technical data **Sedna**

| | |
|------------------|-------------------|
| ■ Total width | 530 mm |
| ■ Total length | 960 mm |
| ■ Total height | 840 mm (870 mm**) |
| ■ Seat height | 540 mm |
| ■ Turning radius | 820 mm |
| ■ Batteries | 2x12V - 11,5Ah |

| | |
|--------------------------|-------------------|
| ■ Driving range | ca. 15 km |
| ■ Power | 270 W |
| ■ Speed | 8 km/h |
| ■ Weight | 27,3 kg (31,61**) |
| ■ Users weight | 115 kg |
| ■ Additional information | CE |



* Driving range under optimal environmental conditions and battery condition, determined according to ISO 7176-4.
 ** Only for the black, Premium model. Measurements in mm, kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible



VERMEIREN India Rehab Private Ltd.

400, Red Sanders, Boulevard, Sri City (SEZ),
 Satyavedu (M), Chittoor District,
 Andhra Pradesh - 517646, India

Tel : +91 8576 666600

www.vermeiren.in
info@vermeiren.in

For more information, contact your local dealer



VERMEIREN

Opening up new horizons in mobility

Subject to technical alterations