

# Jazz S50

## Standard wheelchair

The wheelchair **Jazz S50** is an excellent combination between technicality, comfort and **VERMEIREN** quality. Considered as a standard wheelchair but with lots of benefits.

Ergonomic handgrips, easy to clean nylon upholstery, adjustable arm- and legrests and detachable rear wheels.



## Options

The **Jazz S50** can be equipped as transit model (T30) or extended with several options such as travelwheels, anti-tipping or a B05 armrest. Please consult your dealer or the orderform for a complete and actual list.



**B78 - Anti-tipping wheels**



**B86 - Travelwheels**



**B05 - Adjustable armrest**



**B69 - Foldable backrest**

## Colour(s)



## Technical data

■ Total width	565 / 595 / 615 / 635 / 655 / 675 mm
■ Folded width	290 mm
■ Seat with	390 / 420 / 440 / 460 / 480 / 500 mm
■ Height armrests	220 mm
■ Seat height	500 mm
■ Height backrest	420 mm

## Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

## Quality guaranteed

Regarding quality, the **Jazz S50** complies with the highest norms, laid down by certified test institutes in most European countries.



**T30 - Transit version**

■ Total height	935 mm
■ Seat depth	450 mm
■ Total length	1120 mm
■ Weight	16,6 kg
■ Users weight	130 kg
■ Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

## VERMEIREN India Rehab Private Ltd.

# 400, Red Sanders, Boulevard, Sri City (SEZ),  
Satyavedu (M), Chittoor District,  
Andhra Pradesh - 517646, India

[www.vermeiren.in](http://www.vermeiren.in)  
[info@vermeiren.in](mailto:info@vermeiren.in)



For more information, contact your local dealer

**VERMEIREN**  
we care for you