

Lyna 2

2 in 1 Walker



Firm hand grips



Button press adjustment

With the LYNA 2 you actually get 2 models in 1...

With just one click, it turns from a fixed walker to an articulating walker. The legs are adjustable in height and it has a convenient height to be positioned just above the toilet.

total width	540 mm
total length	470 mm
total height	780 - 960 mm
weight	265 kg
max. users weight	130 kg
colour(s)	inox

Measurement tolerances: ± 15 mm / 1,5kg.
We reserve the right to introduce technical changes.



optional wheels



Maureen



Amy



April



VERMEIREN India Rehab Private Ltd. - # 400, Red Sanders, Boulevard, Sri City (SEZ) - Satyavedu (M), Chittoor District - Andhra Pradesh - 517646, India - www.vermeiren.in - info@vermeiren.in

Model shown only for illustration purposes. Certain features may be optional.

Opening up new horizons in mobility