

# VMS Gynaec Examination Table



VEX-1004

- Strong and Sturdy construction made from mild steel ERW rectangular tubes and machine pressed mild steel CRCA sheet.
- Fitted with rexene covered plain cushioned top with removable U cut.
- Stainless Steel bowl, a pair of adjustable SS lithotomy rods with rexene ankle straps.
- Upper section of the box construction with three separate drawers.
- Height adjustable B.P. apparatus tray provided near head side.
- Built-in sliding step stool.
- Pre-treated and epoxy powder coated finish.



Imaging  
Healthcare  
Signage