

# VMS Manual ICU Bed

Intensive  
Care Unit



The VMS Manual beds combine the full features of height and back adjustment, and TR / reverse TR as well as upper leg adjustment on four separate crank systems. Designed to deliver the VMS quality and safety at an affordable price, these beds are ideally suited for clinics, nursing homes, and hospitals.

V1B-6001



 **VMS** | *Imaging  
Healthcare  
Signage*

# VMS Manual ICU Bed



VIB-6002

- Removable polymer moulded / SS head & foot boards with laminated panels
- Backrest and height adjustment, Upper leg, TR and RTR adjustment on crank mechanism with separate screw mechanism provided at leg end with stainless steel handles
- Four section, perforated mild steel CRCA sheet top
- Mounted on four synthetic 125mm body castors with brakes on two wheels for additional safety
- Rubber buffers on four corners
- Telescopic stainless steel I.V. poles with 4 locations at 4 corners of the bed
- Pre-treated and epoxy powder coated finish.

