

VMS Obstetric Labour Table




VLT-1001



- Stainless steel table top in three sections.
- Stainless steel tubular framework.
- Backrest adjustable on ratchet mechanism.
- Middle Section with "U" cut.
- Leg section can be slid underneath main frame for lithotomy positions.
- Trendelenburg position on hand operated gas spring mechanism.
- Removable stainless steel head and side railings, stainless steel hand grips and removable bowl.
- A pair of rubber padded stainless steel lithotomy rods with rexene ankle straps.
- Supplied in semi knocked down condition.



VMS Obstetric Labour Table



- Handle wound Trendelenburg position : 0 - 10
- Backrest Angle: 0 - 65°
- Height Adjustment: Min: 550mm - Max: 885mm