

## PLAIN BEDS



**VPB-1003 Plain Bed With Removable MS head & foot boards**

Overall approx dimension:  
1980mm (L) x 910mm (W) x 600mm (H). Bed frame approx. dimension: 1980mm (L) x 910mm (W). Uniformly perforated mild steel CRCA sheet top. Framework made from mild steel ERW rectangular tubes. Bed fitted with removable mild steel head & foot boards. "H" type legs fitted with rubber shoes. Pre-treated and powder coated finish. Supplied in knocked down condition.



**VPB-1004 Plain Bed With Un-equal MS head & foot boards**

Overall approx dimension:  
2050mm (L) x 910mm (W) x 600mm (H). Bed frame approx. dimension: 1980mm (L) x 910mm (W). Uniformly perforated mild steel CRCA sheet top. Framework made from mild steel CRCA sheet. Fitted with un-equal height mild steel bows. Tubular legs fitted with rubber shoes. Pre-treated and powder coated finish. Supplied in knocked down condition.



**VPB-1006 Attendant Bed With "H" Type Legs  
(Approx Dimension : 6' x 3')**

Overall approx dimension:  
1830mm (L) x 910mm (W) x 400mm (H). Uniformly perforated mild steel CRCA sheet top. With "H" type legs fitted with rubber shoes. A pair of mattress arrestor. Pre-treated and powder coated finish. Supplied in knocked down condition.

### Optional Accessories for Plain / Attendant Beds:

- VMT-1001 PU Foam Mattress – Single Section Suitable For Plain Bed
- VMT-1003 PU Foam Mattress – Single Section Suitable For Attendant Bed

### Alternate:

**VPB-1007 Attendant Bed With "H" Type Legs  
(Approx Dimension: 6' X 2')**

Overall approx dimension: 1830mm (L) x 610mm (W) X 400mm (H)