

Intensive
Care Unit

VMS Recovery Bed (4-Section)



VRB-5001

Ideally suited for post-operative care, the VMS Four Section Recovery Bed is a dynamic product having the perfect functionalities for stabilising the patient's health condition. Therefore recommended for post-operative recovery rooms and emergency wards as an alternative option.



VMS Recovery Bed (4-Section)



VRB-5002

- Removable polymer moulded / SS head & foot boards with laminated panels
- Backrest and height adjustment, Upper leg, TR and RTR adjustment on crank mechanism with separate screw mechanism provided at leg end with stainless steel handles
- Four section, perforated mild steel CRCA sheet top
- Mounted on four synthetic 125mm body castors with brakes on two wheels for additional safety
- Rubber buffers on four corners
- Telescopic stainless steel I.V. poles with 4 locations at 4 corners of the bed
- Pre-treated and epoxy powder coated finish.



Intensive
Care Unit

VMS Recovery Bed (2-Section)

- Two section, perforated mild steel CRCA sheet top
- Backrest and Trendelenburg and Reverse Trendelenburg adjustment on crank mechanism with separate screw mechanism provided at leg end with stainless steel handles



VRB-5003



VRB-5004

